



Technical Guidance for Laboratories Accredited for EPA Method 200.7

Purpose of This Document

This document provides technical clarification to laboratories performing trace metals analysis by EPA Method 200.7. Some of the critical detail in this method is difficult to interpret.

This document provides guidance to help laboratories implement the necessary steps in the correct order, to ensure the method is followed as intended.

Publication information

This report is available on the Department of Ecology's website at <https://fortress.wa.gov/ecy/publications/SummaryPages/1603028.html>

Ecology's Activity Tracker Code for this study is 10-209.

Website

Laboratory Accreditation: www.ecy.wa.gov/programs/eap/labs

Contact information

Author: Rebecca Wood

Email: rewo461@ecy.wa.gov

Lab Accreditation Unit

Washington State Department of Ecology

PO Box 488, Manchester, WA 98353-0488

7411 Beach Drive East, Port Orchard, WA 98366

Phone: (360) 871-8840

Communications Consultant, Olympia

Phone: (360) 407-6764

Washington State Department of Ecology - www.ecy.wa.gov/

- Headquarters, Olympia (360) 407-6000
- Northwest Regional Office, Bellevue (425) 649-7000
- Southwest Regional Office, Olympia (360) 407-6300
- Central Regional Office, Union Gap (509) 575-2490
- Eastern Regional Office, Spokane (509) 329-3400

Accommodation Requests: To request ADA accommodation including materials in a format for the visually impaired, call Ecology at 360-407-6764. Persons with impaired hearing may call Washington Relay Service at 711. Persons with speech disability may call TTY at 877-833-6341.

Technical Guidance for Labs Accredited for EPA Method 200.7

1. Select suitable wavelengths for analytes of interest. Make sure at least one suitable wavelength has been selected for each analyte of interest. Analyzing two wavelengths is recommended, where possible. The two wavelengths should show similar results. If there is great disparity between wavelengths for the same analyte, this indicates interferences may be affecting the results. Investigate prior to reporting data. It is OK to report from an alternate wavelength as long as supporting Quality Assurance (QA) information is available such as Linear Dynamic Range (LDR), method detection limit, and quality control.
2. Optimize instrument and set background points (Method 200.7 section 4.1.1 and Method 200.5 section 4.1). Background points should be set prior to determining LDR and Inter-element Corrections (IECs). Avoid locating background points within regions that are impacted by spectral interferences.
3. Identify interferences, and ensure that they are accounted for in the IEC correction routine.
4. Perform LDR study. LDR must be within 10% of the true value. LDRs are part of the Initial Demonstration of Performance (IDP), which is mandatory (Method 200.7 and Method 200.5 section 9.2.2).
5. Establish Spectral Interference Correction (SIC) Routine (Method 200.7 section 10.4 and Method 200.5 section 4.1 note). Run single-element solutions which must be within the LDR. Typically 50 to 100 ppm (Method 200.7 section 4.1.4). Single element solutions should be analyzed for all analytes of interest and all interferences at the selected wavelengths.
6. Analysis of the single element solutions will generate apparent false positive or false negative peaks. These will be used to calculate K factors which will be stored in the instrument software. Divide the apparent concentration by the measured concentration of the interference to calculate the k value. Some software is programmed to perform these calculations and store the data.
7. Once established, the IEC routine must be verified. A daily SIC solution should be analyzed including all major cations present in the samples. Typically Al, Mg, Ca, and Fe for most ambient, drinking waters and ground waters. In some samples such as influents and waste samples, other interferences may need to be included in the SIC solution. Daily SIC solutions should include high concentrations such as 50 to 100 ppm.
8. Acceptance criteria for daily SIC may be set at +/- LLOQ/PQL.
9. The complete IEC correction routine should be verified annually.
10. Perform Quality Control Sample (QCS) study (Method 200.7 section 9.2.3 and Method 200.5 section 9.2.4). Second source standard +/- 5% of true value.

11. Perform MDL studies annually (Method 200.7 section 9.2.4 and Method 200.5 section 9.2.5). MDL studies should be wavelength specific and be taken through the entire preparation and analysis.

12. Ensure calibration standards are prepared as described in Table 3 of Method 200.7. When all elements are added to one standard, spectral interference error may occur in the calibration curve, since the correction routine is not applied until after calibration.

The use of internal standards is recommended for axial view, especially for water samples that are not particularly clean, or soils.